

Sefar Solutions for Woundcare

SEFAR MEDIFAB fabric as a single-layer overlay or SEFAR MEDITEX ACC as a multi-layer overlay allow unsurpassable professional treatment of different wounds.



Product Features

Function

SEFAR MEDIFAB and SEFAR MEDITEX ACC allow the treatment of different types of wounds. Whether as a contact layer for low-pressure wound therapy, as a supporting fabric in skin grafts or as a carrier for drugs and ointments, these non-adherent products are characterized by having well-defined surface properties and a homogeneous structure.

By choosing a specific yarn or equivalent coating, numerous other functions (e.g., antibacterial) can be achieved, which facilitates the development of new high performance wound dressings.

DOWNLOADS

Brochure: Medical Components (PDF 3997 kb)

Customer Info: SEFAR MEDIFAB Fabrics are Low Protein Binding (PDF 208 kb)

Your Benefits

- Non-stick surfaces
- Defined pore/mesh size
- High throughput rates
- Non-shedding due to monofilament yarn
- Unsurpassed, two-dimensional stability
- Can be combined with a wide range of absorbers and nonwovens
- Can be manufactured to customers' specific requirements

Sefar AG

Hinterbissastrasse 12
 9410 Heiden
 Switzerland

Phone +41 71 898 5700
 Fax +41 71 898 5721


info@sefar.com

[Go to product page](#)

Product range



■ SEFAR MEDIFAB



■ SEFAR MEDITEX
ACC

For further information and fabric samples, please select the preferred contact option on the right.

Locations



Sefar AG

Töberstrasse 4
9425 Thal – Switzerland
Phone: +41 71 898 5700
Fax: +41 71 886 3504

■ E-Mail



Sefar AG

Hinterbissastrasse 12
9410 Heiden
Switzerland
Phone: +41 71 898 5700
Fax: +41 71 898 5721

■ E-Mail

Sefar AG

Hinterbissastrasse 12
9410 Heiden
Switzerland

Phone +41 71 898 5700
Fax +41 71 898 5721

info@sefar.com

[Go to product page](#)